



- 1** Cage type II
- 2** Filter hood for cage type II
- 3** Stainless steel cover for filter hood, inside-located
- 4** Stainless steel cover with falling wire clip
- 5** Stainless steel cover with clamp closing

Cage type II

The cages are available in different materials depending on the application. The corresponding stainless steel covers are provided with a feeding place across the entire width for holding the watering bottle and the food. The feeding place can be separated by a divider of stainless steel. On application of the filter hood, the inside-located stainless steel cover with the integrated divider must be used. The filter hood is provided with a filter sheet of class G4 and can be delivered optionally with fixing clamps and a peripheral silicone gasket.

As further options you will find within this group different feeders, watering bottles, drinking tops and label holders.

Type	Material	Dimensions W x H x D	Floor surface	Weight	Code-no.
Cage size II	Polycarbonate PC 3108	207 x 140 x 265 mm	363 cm ²	approx. 440 g	II-1264.PC
Cage size II	Polysulphone PSU	207 x 140 x 265 mm	363 cm ²	approx. 440 g	II-1264.PSU
Cage size II	Polypropylene PP	207 x 140 x 265 mm	363 cm ²		II-1264.PPN
Cover with falling wire clip	St. steel AISI 304	Grid distance 7 mm			II-1264.F
Cover with clamp closing	St. steel AISI 304	Grid distance 7 mm			II-1264.K
Divider	St. steel AISI 304				1264-TR
Cover for filter hood, inside-located with divider	St. steel AISI 304	Grid distance 7 mm			II-1264.FIDE
Filter hood	Polycarbonate PC 3108				II-1264.FIHA
Filter hood	Polysulphone PSU				II-1264.FIHA
Floor	St. steel AISI 304				II-1264-LB