



- 1** Bottle 900 ml
- 2** Bottle 700 ml
- 3** Bottle 500 ml
- 4** Bottle 300 ml
- 5** Bottle 150 ml
- 6** Bottle 250 ml with groove
- 7** Bottle 400 ml with groove
- 8** Bottle 700 ml with groove

Watering bottles

The Ehret watering bottles are made of high quality polycarbonate, some also of polysulphone. Round as well as square models are available for selection. All the round bottles have a clearly visible scaling. In case the bottles shall be autoclaved permanently in filled condition and including the drinking top, the models with groove and silicone ring must be used.

All the bottles that are named here can be combined with the Ehret drinking tops.

Type	Bottom Dimensions	Shape	Material	Code-no.
150 ml	Ø 54 mm	Round	PC	FL-150-54
250 ml	Ø 54 mm	Round	PC	FL-250-54
300 ml	Ø 84 mm	Round	PC	FL-300-84
500 ml	Ø 84 mm	Round	PC	FL-500-84
700 ml	Ø 84 mm	Round	PC	FL-700-84
900 ml	Ø 84 mm	Round	PC	FL-900-84
250 ml / with groove	Ø 54 mm	Round	PC	FL-250-NUT
400 ml / with groove	70 x 70 mm	Square	PC	FL-400-NUT
530 ml / with groove	70 x 70 mm	Square	PC	FL-530-NUT
700 ml / with groove	70 x 70 mm	Square	PC	FL-700-NUT
300 ml / with groove	Ø 58 mm	Round	PSU	FL-300-NUT- Ø58-PSU
500 ml / with groove	Ø 66 mm	Round	PSU	FL-500-NUT- Ø66-PSU
400 ml / with groove	70 x 70 mm	Square	PSU	FL-400-NUT-quadr.-PSU
530 ml / with groove	70 x 70 mm	Square	PSU	FL-530-NUT-quadr.-PSU
750 ml / with groove	73 x 73 mm	Square	PSU	FL-750-NUT-quadr.-PSU