



- 1** Stretched nozzle
- 2** Nozzle 10 mm
- 3** Nozzle 25 mm
- 4** Nozzle 40 mm
- 5** Ball nozzle 40 mm
- 6** Ball nozzle 100 mm, straight
- 7** Nozzle 100 mm, curved
- 8** Ball nozzle 140 mm, curved

Drinking tops

The water supply for small rodents, rabbits as well as for primates can be covered with the available drinking tops. The drinking tops can be used with or without silicone ring.

All drinking tops are made of high quality stainless steel AISI 304 and are therefore autoclavable.

| Opening diameter | Description | Code-no. |
|------------------|------------------------------|------------|
| 0,7 mm | Stretched nozzle | 2004-96 |
| 1,5 mm | Nozzle 10 mm | 97-10-7 |
| 1,5 mm | Nozzle 25 mm | 97-25-7 |
| 1,5 mm | Nozzle 40 mm | 97-40-7 |
| 5,4 mm | Ball nozzle 40 mm | 97-40-K8 |
| 5,4 mm | Ball nozzle 100 mm, straight | 97-100-K8 |
| 5,4 mm | Ball nozzle 100 mm, curved | 97-100-GK8 |
| 5,4 mm | Ball nozzle 140 mm, curved | 97-140-GK8 |