



AIR PURIFIER QRJ – 128

MOBILE OR FIXED

The QRJ-128 air purifier offers protection with horizontal laminar flow. It is an ideal system to protect from contamination work in the room where it is required. This unit is widely used in hospitals, research laboratories, clean rooms, etc.

CHARACTERISTICS

- Offers air rating 100 (ISO 5) at the outlet
- The unit is harmless, innocuous, low noise, and high efficiency performing the most effective, fast, and complete air purification.
- Absorbs all aerosols and contaminants from the room or laboratory.
- Reduces room contamination through air recirculation.
- Output flow regulation.
- High efficiency HEPA filter for 0.3 micron particles.
- Its special design allows its installation in multiple positions. It can be mobile or fixed to wall or ceiling.
- Low noise level < 58 dB



<< Data display
Remote control >>



PRINCIPAL OFFICE

José Echegaray, 8 Local 1-3
28100 – Alcobendas – Madrid - Spain
Tel. +34 916 624 473 Fax.+34 916 624 469
www.tdi.es ventas@tdi.es

DELEGATIONS

Catalunya. Tel. +34 676 351 264
Galicia. Tel. +34 639 486 955
Andalucía. Tel. +34 639 481 756

AIR PURIFIER QRJ – 128

MOBILE OR FIXED

CHARACTERÍSTICS (cont)

- Includes internal UV decontamination
- Programming capacity of the start-up and shutdown of the flow system, in up to 3 separate cycles
- UV decontamination programming capability in up to three separate cycles
- Possibility of mixing programming cycles. For example.

9:00 Filtering ignition

14:00 Uv decontamination ignition

16:00 Turning off the leak

18:00 off of UV decontamination

Etc...

- Remote control and programming via IR remote control, for units located in low accessibility locations (ceilings, high walls, etc...)
- Possibilities of work on mobile support (included), Hanging from the ceiling (single or forming joint systems) and recessed in the wall.

TECHNICAL SPECIFICATIONS

MODEL	QRJ-128
Dimensions	900 x 120 x 1500 mm
Filtration surface	720 x 570 mm
Air Velocity	>0.25 m/s
Weight	58 Kg
Airflow	350 m ³ /h
Noise level	< 58 db
Power	220 v , 50 Hz
Air Clasificación	Class 100

ORDER INFORMATION

Ref	Description
425-QRJ-128	Air Purificator / Mobile Operating Unit horizontal flow QRJ-128

PRINCIPAL OFFICE

José Echegaray, 8 Local 1-3
28100 – Alcobendas – Madrid - Spain
Tel. +34 916 624 473 Fax.+34 916 624 469
www.tdi.es ventas@tdi.es

DELEGATIONS

Catalunya. Tel. +34 676 351 264
Galicia. Tel. +34 639 486 955
Andalucía. Tel. +34 639 481 756