

# CONSTANT TEMPERATURE INCUBATOR

## Characteristics

- LCD Display.
- Microcomputer temperature controller
- Double door design and magnetic door sealing
- The wheels are equipped for easy movement.
- The working chamber is designed with rounded corners that allows easy to clean.
- Protection against elevated temperatures with temperature sensor malfunction alarm.
- Parameters can be saved automatically when it shuts down or a system error occurs.



**425-ITC-2P-50/80/160/270 II**



**425-ITC-2P-35L**

### OFICINA PRINCIPAL

### DELEGACIONES

## Specifications

MODEL	ITC-2P-35L	ITC-2P-50L	ITC-2P-80L	ITC-2P-160L	ITC-2P-270L
Internal size (width x depth x heigth) (mm)	340 x 320 x 320	350 x 348 x 400	410 x 405 x 510	505 x 455 x 700	600 x 600 x 750
External size (mm)	620 x 490 x 490	510 x 548 x 740	570 x 590 x 850	673 x 645 x 1033	768 x 795 x 1083
Package size (mm)	740 x 600 x 670	590 x 650 x 930	730 x 700 x 1030	800 x 800 x 1210	930 x 940 x 1250
Weight (kg)	31	50	75	100	150
Capacity (L)	35L	50L	80L	160L	270L
Temperature range	RT + 5 - 65°C	RT + 5 - 60°C			
Temperature Fluctuation	± 0,3°C	± 0,3°C			
Time range	0-5999 min	1-9999 min			
Shelves	2 adjustable			3 adjustable	
Power supply	220 V/50 Hz				
Consumption (W)	200	437	502	580	688

## Order Information

Ref	Description
425-ITC-2P-35L	CONSTANT-TEMPERATURE INCUBATOR 2 DOORS 35L
425-ITC-2P-50L	CONSTANT-TEMPERATURE INCUBATOR 2 DOORS 50L
425-ITC-2P-80L	CONSTANT-TEMPERATURE INCUBATOR 2 DOORS 80L
425-ITC-2P-160L	CONSTANT-TEMPERATURE INCUBATOR 2 DOORS 160L
425-ITC-2P-270L	CONSTANT-TEMPERATURE INCUBATOR 2 DOORS 270L