

→ Universal racks



Laboratory Animal Technology: Cages and accessories for rodents





1



2

1 Universal rack THF 3057

2 Shelf made of round tubes

Universal racks

The sturdy universal racks consist of modules, which can be combined. They are made of stainless steel AISI 304 and are completely autoclavable. The glass bead blasted surface makes cleaning and disinfecting easy.

The following options are available:

- shelves made of round tube material, Ø 6 mm
- shelves made of metal sheet
- welded shelves
- fitted with feet
- fitted with autoclavable wheels
- screw-connected shelves

The universal racks can be delivered in various dimensions and as corner-fitting combination.

Type	Material	Dimensions W x D x H	Wheels	Design	Code-no.
Universal rack	St. steel AISI 304	1280 x 600 x 1790 mm	4, 2 of them with brakes	Peripheral angle material with shelves of round tubes, Ø 6 mm	THF 3057
Universal rack	St. steel AISI 304	1280 x 600 x 1790 mm	4, 2 of them with brakes	Peripheral angle material with shelves of round tubes, Ø 6 mm	THF 3057
Universal racks for combination	Please let us know your requirements. We will gladly advise you.				THF 3057