




→ Racks for plastic cages



**Laboratory Animal Technology:
Cages and accessories for rodents**



1 Rack for plastic cages
Code-no. 2152/E, hairpin

2 Rack for plastic cages
Code-no. 2154/E

3 Rack for plastic cages
Code-no. 2494/E

4 Detail:
hairpin bracket

5 Detail:
U-profile bracket

Racks for plastic cages

The racks for plastic cages are available for the prevalent cage types I L to IV. The cleaning and disinfecting of the stainless steel AISI 304 racks is facilitated by the glass bead blasted surface. Optionally hairpin brackets as well as U-profile brackets are available. The 4 wheels, Ø 100 mm, are autoclavable, 2 of them are fitted with brakes.

Options:

- wheel Ø 125 mm
- bumpers
- cover against incidence of light

Rack for:	Construction	Dimensions W x H x D	Code-no. Stainless steel AISI 304 with hairpin bracket	Code-no. st. Stainless steel AISI 304 with U-profile bracket
Cages type I L	Single-sided	1296 x 1765 x 550 mm	THF 3328	THF 3321-2/E
Cages type I L	Double-sided	1296 x 1765 x 710 mm	THF 3329	THF 3322-2/E
30 cages type II	Single-sided	1296 x 1765 x 550 mm	THF 2364	THF 2151/E
35 cages type II	Single-sided	1296 x 1765 x 550 mm	THF 2365	THF 2152/E
60 cages type II	Double-sided	1296 x 1765 x 710 mm	THF 2366	THF 2153/E
70 cages type II	Double-sided	1296 x 1765 x 710 mm	THF 2367	THF 2154/E
30 cages type II L	Single-sided	1296 x 1765 x 550 mm	THF 3145-5	THF 3101-5(3)/E
35 cages type II L	Single-sided	1296 x 1765 x 550 mm	THF 3145-6	THF 3101-6(3)/E
60 cages type II L	Double-sided	1296 x 1765 x 710 mm	THF 3145-7	THF 3101-7(3)/E
70 cages type II L	Double-sided	1296 x 1765 x 710 mm	THF 3145-8	THF 3101-8(3)/E
20 cages type III	Single-sided	1296 x 1765 x 550 mm	THF 2796	THF 2580/E
24 cages type III	Single-sided	1296 x 1765 x 550 mm	THF 2368	THF 2155/E
40 cages type III	Double-sided	1296 x 1765 x 780 mm	THF 2797	THF 2795/E
48 cages type III	Double-sided	1296 x 1765 x 780 mm	THF 2369	THF 2156/E
15 cages type IV	Single-sided	1294 x 1780 x 590 mm	-	THF 2494/15/E
18 cages type IV	Single-sided	1294 x 1780 x 590 mm	-	THF 2494/18/E
20 cages type IV	Single-sided	1717 x 1780 x 590 mm	-	THF 2494/20/E
24 cages type IV	Single-sided	1717 x 1780 x 590 mm	-	THF 2494/24/E