



Labor- und Pharmatechnik

→ **Cage type IV**



**Laboratory Animal Technology:  
Cages and accessories for rodents**




EHRET LIFE SCIENCE SOLUTIONS




Labor- und Pharmatechnik




1



2



3



4



5

1 Cage type IV

2 Stainless steel cover without feeding place

3 Floor

4 Stainless steel cover, raised design

5 Stainless steel cover with clamp closing

## Cage type IV

The cages are available in different materials depending on the application. The corresponding stainless steel covers are provided with a feeding place across the entire width for holding the watering bottle and the food. The feeding place can be separated by a divider of stainless steel. The height-increased stainless steel cover is used for the keeping of rats – thereby an optimal cage height is achieved and the animals have the possibility to communicate between the cages. Furthermore, stainless steel covers for the husbandry of guinea pigs are available.

As further options you will find within this group different feeders, watering bottles, drinking tops and label holders.

Type	Material	Dimensions W x H x D	Floor surface	Weight	Code-no.
Cage size IV	Polycarbonate PC 3108	380 x 200 x 590 mm	1815 cm <sup>2</sup>	approx. 2100 g	40000-PC
Cage size IV	Polysulphone PSU	380 x 200 x 590 mm	1815 cm <sup>2</sup>	approx. 2100 g	40001-PSU
Cover with falling wire clip	St. steel AISI 304	Grid distance 7 mm			40050
Cover with clamp closing	St. steel AISI 304	Grid distance 7 mm			40055
Cover with clamp closing, raised design	St. steel AISI 304	Grid distance 7 mm + 50 mm height			40051
Cover for guinea pigs	St. steel AISI 304	Grid distance 16 mm			40070
Divider	St. steel AISI 304				40060
Floor	St. steel AISI 304				40065