




Labor- und Pharmatechnik

→ Bio. A.S. Cage type II L



Laboratory Animal Technology: Air Flow Systems



EHRET LIFE SCIENCE SOLUTIONS



Labor- und Pharmatechnik



1 Cage type II L 2 Stainless steel cover for filter hood, inside-located 3 Filter hood, bottle outside 4 Filter hood, bottle inside

Cage type II L

Bio. A.S. cage type III L is optimised to a maximum of tightness and is provided with a peripheral silicone gasket. Contamination of the rack and the air conduction pipes is avoided by the ventilation cover designed by EHRET. Furthermore a homogenous air distribution is achieved at low flow velocities. The system is available with drinking bottle outside as well as inside. The filter hoods can optionally be equipped with a CO₂ membrane. (see image on Bio. A.S. cage type IL)

EHRET Bio. A.S. cages guarantee maximum biological safety, easy and ergonomic handling and conformity with the legal requirements.

Type	Material	Dimensions W x H x D	Code-no.
Cage type II L	Polycarbonate PC Polysulphone PSU Polyetherimid PEI	207 x 140 x 365 mm floor space 530 cm ²	20000-PC 20001-PSU 20002-PEI
Stainless steel cover, inboard, with collapsible dividing sheet	Stainless steel AISI 304 electropolished	Grid distance 7 mm	20550I
Stainless steel cover, inboard, with divided food place	Stainless steel AISI 304 electropolished	Grid distance 7 mm	20551
Bio. A.S. filterhood type II Long, drinking bottle outside, complete standard version consisting of: filter hood, ventilation cover, silicone gasket and fastening clamps	Polycarbonate PC Polysulphone PSU Polyetherimid PEI	Height 45 mm	20500-PC 20501-PSU 20502-PEI
Bio. A.S. filterhood type I Long, drinking bottle outside, complete standard version as described above, additionally equipped with: closures in supply / exhaust air	Polycarbonate PC Polysulphone PSU Polyetherimid PEI	Height 45 mm	20520-PC 20521-PSU 20522-PEI
Bio. A.S. filterhood type I Long, drinking bottle outside, complete standard version as described above, additionally equipped with: closures in supply / exhaust air and CO ₂ membrane PALLFLEX TL 1865 (THF3302-2) with covering	Polycarbonate PC Polysulphone PSU Polyetherimid PEI	Height 45 mm	20530-PC 20531-PSU 20532-PEI
Options:			
Car holder DIN A7 with short clip	Stainless steel AISI 304		THF1562-2
Card holder DIN A7 with T-clip, portrait format, collapsible	Stainless steel AISI 304		THF1562-8
Card holder DIN A7 with T-clip, landscape format, collapsible	Stainless steel AISI 304		THF1562-9
Card holder DIN A7 with T-clip, portrait format, to be attached to the filterhood	Stainless steel AISI 304		THF1562-16

1* consisting of 1 filter hood, 1 set ventilation cover, 1 set filter fleece strips, 1 silicone gasket, 2 fixing clamps